CONTENT

EN (English) Instruction Manual Zfx™ Screws (incl. GenTek™ Screws)	2
BG (Български) Ръководство за употреба Винтове Zfx™ (включително винтове GenTek™)	4
CS (Česky) Návod k použití Šrouby Zfx™ (včetně šroubů GenTek™)	6
DA (DANISH) Brugervejledning Zfx™-skruer (inkl. GenTek™-skruer)	8
DE (Deutsch) Gebrauchsanweisung Zfx™ Schrauben (inkl. GenTek™ Schrauben)	10
EL (Ελληνικά) Εγχειρίδιο οδηγιών Βίδες Zfx™ (συμπεριλ. βίδες GenTek™)	12
ES (Español) Manual de instrucciones Tornillos Zfx™ (incl. tornillos GenTek™)	14
ET (Eesti) Kasutusjuhend Zfx™ Kruvid (sh GenTek™-i kruvid)	16
FI (Suomi) Käyttöohje Zfx™-ruuvit (ml. GenTek™-ruuvit)	18
FR (Français) Mode d'emploi Vis Zfx™ (Inclut les vis GenTek™)	20
HR (Hrvatski) Priručnik Vijci Zfx™ (uklj. vijke GenTek™)	22
HU (Magyar) Kezelési útmutató Zfx™ Csavarok (GenTek™ csavarok is beleértve)	24
IT (Italiano) Manuale di istruzioni Viti Zfx™ (comprese viti GenTek™)	26
LT (Lietuvių k) Naudojimo instrukcija „Zfx™“ varžtai (įskaitant „GenTek™“ varžtus)	28
LV (Latviski) Lietošanas rokasgrāmata Zfx™ skrūves (kā arī GenTek™ skrūves)	30
NL (Nederlands) Gebruikershandleiding Zfx™-schroeven (incl. GenTek™-schroeven)	32
NO (Norsk) Instruksjonshåndbok Zfx™ skruer (inkl. GenTek™ skruer)	34
PL (Polski) Instrukcja stosowania Śruby Zfx™ (łącznie ze śrubami GenTek™)	36
PT (Português) Manual de Instruções Parafusos Zfx™ (incl. parafusos GenTek™)	38
RO (Română) Manual cu instrucțiuni Șuruburi Zfx™ (inclusiv șuruburi GenTek™)	40
SK (Slovak/Slovakian) Návod na použitie Skrutky Zfx™ (vrátane skrutiek GenTek™)	42
SL (Slovenščina) Navodila za uporabo Vijaki Zfx™ (vklj. z vijaki GenTek™)	44
SV (Svenska) Bruksanvisning Zfx™ Skruvar (inkl. GenTek™-skruvar)	46

Indications


Intended Use / Intended Purpose

The Zfx™ Screw and GenTek™ Screw are intended to screw Ti-Base and direct screwed prosthetic components on cement-retained single-unit, multi-unit, and full-arch prosthetic restorations over implant or analog. The Zfx™ Screw and GenTek™ Screw is not intended for restorations requiring an angle more than 20° to the axis of the implant. Please also considerate the IFUs for our Zfx™ Titanium Base / Zfx™ Ti-Abutment Blank and GenTek™ Titanium Base / GenTek™ Titanium pre-milled Abutment Blank, located on the Zfx website.


Device description

Zfx™ Screws and GenTek™ Screws are manufactured from biocompatible titanium alloy grade 5 ELI (TiAl6V4 ELI). Zfx™ Screws and GenTek™ prosthetic screws are for single use only. GenTek™ Gold-Tite Screws have a surface in pure Gold-plating. Zfx™ Screws and GenTek™ prosthetic screws are provided non-sterile and intended to be cleaned and sterilized prior to placement into the patient's mouth. A new screw should be used when assembling a prosthesis for the first time and for every check thereafter.

Screws are available with two different screw driver engagements:

Hexalobular Drive

Hex Drive

Please make sure to use the correct corresponding Screwdriver Tip (Hex or Hexalobular).

GenTek™ Screws

The GenTek™ Screws are available for genuine connections for the following Zimmer Biomet Dental implant systems: Certain®, External Hex, TSV®/Trabecular Metal™, and Eztetic®. For specific dental implant compatibility, please refer to the information on the individual product label. GenTek™ Screws are compatible with Zimmer Biomet Dental prosthetic tooling as outlined in the GenTek™ Components Catalogs – located on the Zfx website. Please make sure to use the correct and corresponding Screwdriver Tip.

GenTek™ Try-In screws

The GenTek™ Try-In screws are available for genuine connections for the following Zimmer Biomet Dental implant systems: Certain®, External Hex and Low Profile. Try-In screws are for exclusively laboratory use and processing of the abutment.

Zfx™ Screws

The Zfx™ Screws are available for several systems. Please refer to the Zfx Components Catalog located on the Zfx website.

Screwdriver

The correct Screwdriver Tip must be used (see references to document on website) to avoid damage at the screw head. Make sure that the Tip fits correctly with the corresponding torque wrench.

Recommended Torque

For each screw type a recommended torque value for the final placement is defined. The recommended torque value regarding the used screw must be applied to assure a correct fitting of the TiBase or Abutment on the implant. Don't exceed the applied torque over this value to avoid damaging the TiBase/ Abutment and the implant! For the recommended torque value and further information please refer as well to the relative IFUs of the TiBase, Abutment Blank chosen for the project.

Cleaning and sterilization

The Zfx™ Screws and GenTek™ Screws are non-sterile. Prior to placement of the Zfx™ and GenTek™ TiBase / Screw restoration and prosthetic screw into the implant, the Zfx™ Screw and GenTek™ TiBase / Screw assembly must first be cleaned and sterilized.

Step 1: Manual Cleaning

1. Rinse the instruments under running cold tap water for 2 minutes.
2. Hand wash instruments in a bath with demineralized water for 5 minutes. Remove debris with a soft brush. Cavities should be rinsed with a syringe. This step should be repeated until all debris is removed.
3. Clean the instrument in an ultrasonic bath using a mildly alkaline cleaner (1.5%) at a frequency of at least 40 kHz for 20 minutes.
4. Rinse thoroughly with demineralized water for 1 min.

Step 2: Sterilization

One of the following sterilization methods is recommended to be used:

- × **Method 1:**
Full cycle pre-vacuum steam sterilization at a temperature of 134 °C (273.2 °F) for 3 minutes, 4 pulse minimum.
- × **Method 2: (alternative to method 1)**
Full cycle gravity steam sterilization at 134 °C (273.2 °F) for 3 minutes minimum.

Step 3: Drying

After sterilization treatment it is recommended to dry the products at a temperature of maximum 134 °C (273.2 °F) for 10 minutes. The products are ready for use 5 minutes after completion of the drying process.

Note: Follow the loading instructions issued by the sterilizer manufacturer.

Known Side Effects

Allergies to the alloy or contents of the alloy may very rarely occur.

Instruction Manual

Zfx™ Screws

(incl. GenTek™ Screws)

Version: 4 / 07.2023
Compiled at: 07.07.2023
Effective date: 07.07.2023
Replaces version: 3 / 01.2022

Warnings

- × The Zfx™ Screw and GenTek™ Screw must only be used and handled by dental professionals who have approved training in dentistry.
- × The Zfx™ Screw and GenTek™ prosthetic screw are small and, therefore, they must be handled with caution to avoid swallowing or inhaling by the patient.
- × The use of the prosthetic screw torque higher than the recommended value may result in damage to the Zfx™ Ti-Base and GenTek™ TiBase, the prosthetic screw and/or the implant. The use of torque values lower than those recommended may result in loosening of the Zfx™ and GenTek™ TiBase and/or the prosthetic restoration.
- × The screwdriver must be in good condition. A worn screwdriver can lead to over-rotation of the hex in the screw head which can result in the screw not being fully tightened and/or not being able to be removed.
- × For any lab work and any try-in of the GenTek™ TiBase, a lab screw must be used. The included prosthetic screw may only be used to permanently fix the finished restoration.

Explanation of symbols

Non sterile

Do not use if package is damaged

Manufacturer Information: Address of legal manufacturer facility

Don't use twice

Recommended torque to be applied for final placement of the included screw

Recommended to fix the Try-In Screw hand-tight

LOT-Number, Symbol followed by Production Date of the lot

Medical Device

Reference Number, Symbol followed by Manufacturer's Item Number

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T +49 (0) 8131 / 33 244 - 0
F +49 (0) 8131 / 33 244 - 10
www.zfx-dental.com

Incident Reporting / Complaints

For Zfx™ products please send to
qm@zfx-dental.com

For GenTek™ products please send to
emeacomplaints@ZimVie.com

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Zfx GmbH