

# Micro 10<sup>+</sup>



**Economical, aldehyde-free concentrated solution for the cleaning and disinfection of medical devices (immersible) and dental instruments before sterilisation, such as root-canal instruments and dental burs.**

## Key product features

- Very economical 2% dilution
- Contact time: 15 minutes
- Suitable for all instruments
- Contains corrosion inhibitors to protect instruments against rusting
- Suitable for ultrasonic cleaners to boost cleaning effect
- VAH/DGHM certificated

# Micro 10<sup>+</sup>


## Protocol for use

### Before use

Wear protective gloves, protective clothing, eyes/face protection.



Dilute Micro 10<sup>+</sup> in a soaking tank at 2% (20ml per 1 litre of water).



Immerse the instruments in the solution and leave for 15 minutes.



Carefully remove the instruments from the solution and inspect. Brush if necessary to remove any remaining visible debris.



Rinse instruments under running water or by immersion in clean water.



Dry instruments with a clean disposable paper towel. Pack the instruments and prepare for sterilisation.


### Recommended dilutions

2% (20ml in 1L water) for all dental instrumentation  
10% (100ml in 1L water) for small instruments,  
e.g. canal instruments, diamond burs, etc.

### Precautions

Dangerous: respect the precautions for use (established according to the European regulations in force for classification and labelling of chemicals). Before each use, read the label and the information concerning the product. Storage between +5°C and +25°C. Reserved for professional use. Medical device CE 0459.

## Efficacy



\* PRV: Surrogate of Hepatitis B

\* BVDV: Surrogate of Hepatitis C

The Robert Koch Institute recommends the use of disinfectants claiming limited virucidal activity to interrupt the infection-chain. According to the RKI, limited virucidal activity (effective against enveloped viruses) corresponds to an activity against Vaccinia virus and BVDV.

### Conforms to European standards

Bactericide: EN 13727, EN 14561

Yeasticide: EN 13624, EN 14562

Active on viruses according to EN 14476

### Product composition

100 g of Micro 10<sup>+</sup> contains 9 g of alkylbenzyltrimethylammonium chloride, amphoteric and non-ionic surfactants, complexant, corrosion inhibitor and additives.

### Physical data

Colour: Pale yellow

Perfume: Without perfume

pH: 7.5 (diluted at 2%)

Shelf Life: Concentrate 3 years

Diluted, unused 28 days

Diluted, used 7 days depending on soil quantity

Container: Made of recyclable HDPE material

### Packaging

1L	2,5L	5L
✓	✓	✓

### Accessories

2L Soaking tank

Soaking tank for burs

