


EN  
RD614

MAVIG




**Note:** In order to meet hygiene requirements, we recommend that you use our washable hygienic cover made from special ComforTex® MF material "Titan" colour.



With interchangeable hook-and-loop fastener



Children's size in Happy design



Available optional hygienic cover

**Models:** RD615E Adult patient thyroid and sternum protection  
RD6145K Children patient thyroid and sternum protection

Made in Germany



## RD614 – Thyroid and Sternum Protection for Patients

The thyroid is a radiation-sensitive organ of the human body and therefore in particular need of protection. The MAVIG thyroid sternum protection for patients is optimally cut and has the best fit, and, of course, it also takes all radiation protection aspects into account.

An uncomplicated magnetic or hook-and-loop fastener, which can also be replaced as an option, provides sufficient scope for adjustment.

- **Protective Material** Standard lead or Leadfree
- **Lead Equivalent** 0.50 mmPb
- **Magnetic/hook-and-loop fastener**
  - Adults Standard size fits neck sizes of 32 - 45 cm
  - Children Standard size fits neck sizes of 28 - 36 cm
- **Interchangeable hook-and-loop fastener**
  - Adults Standard size fits neck sizes of 37 - 50 cm
  - Children Standard size fits neck sizes of 28 - 36 cm
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Plus (designs only), available in the colours:  
*Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or designs: Happy, Nautilus*

**MAVIG GmbH**  
X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com