















**8280-100**  
 to  
**8294-100**
**MEDICOM® DRY-BACK® DENTAL BIBS**
**FEATURES & BENEFITS**

Medicom® Dry-Back® 3-ply dental bibs range is made of 2 layers of embossed cellulose and 1 layer of waterproof polyethylene material. This allows maximum protection with horizontal embossed pattern and superior absorption thanks to additional layer of cellulose.


**PRODUCT DESCRIPTION**

Product name: Dry-Back® Dental Bibs  
 Product type: Single use, non sterile  
 Material: Embossed cellulose, Polyethylene  
 Grammage: 33 g/m<sup>2</sup>  
 Folding: M  
 Size: 33 x 45,5 cm (± 0,5 cm)  
 Unit weight: 6 g (± 10 %)  
 Packaging: Dispenser box of 100 units



Available in 13 colours:	Lavender	Blue	White	Yellow	Green	Dusty Rose	Peach
							
	8280-100	8282-100	8283-100	8284-100	8285-100	8286-100	8288-100
	Aqua	Dark Blue	Orange	Black	Lime	Turquoise	
							
	8289-100	8290-100	8291-100	8292-100	8293-100	8294-100	

Origin: Europe

**INDUSTRIES**

Medical / Dental

**CERTIFICATION**

Manufactured under an ISO 9001 certified quality system.

8280-100  
to  
8294-100

## MEDICOM® DRY-BACK® DENTAL BIBS

### STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C.

Shelf life of the product: 5 years

### LOGISTICAL INFORMATION



#### Outer case specifications

Item	Size mm	Gross weight kg	QTY/pallet	EAN code
8280-100	370 x 390 x 260	4,26	48 cartons	8720249902805
8282-100	370 x 390 x 260	4,26	48 cartons	8720249902829
8283-100	370 x 390 x 260	4,26	48 cartons	8720249902843
8284-100	370 x 390 x 260	4,26	48 cartons	8720249902867
8285-100	370 x 390 x 260	4,26	48 cartons	8720249902881
8286-100	370 x 390 x 260	4,26	48 cartons	8720249902904
8288-100	370 x 390 x 260	4,26	48 cartons	8720249902928
8289-100	370 x 390 x 260	4,26	48 cartons	8720249902942
8290-100	370 x 390 x 260	4,26	48 cartons	8720249902966
8291-100	370 x 390 x 260	4,26	48 cartons	8720249902980
8292-100	370 x 390 x 260	4,26	48 cartons	8720249903000
8293-100	370 x 390 x 260	4,26	48 cartons	8720249903024
8294-100	370 x 390 x 260	4,26	48 cartons	8720249903048



#### Dispenser box specifications

Size mm	QTY/outer case	EAN code
125 x 175 x 240	6 boxes of 100 units	8720249902812
125 x 175 x 240	6 boxes of 100 units	8720249902836
125 x 175 x 240	6 boxes of 100 units	8720249902850
125 x 175 x 240	6 boxes of 100 units	8720249902874
125 x 175 x 240	6 boxes of 100 units	8720249902898
125 x 175 x 240	6 boxes of 100 units	8720249902911
125 x 175 x 240	6 boxes of 100 units	8720249902935
125 x 175 x 240	6 boxes of 100 units	8720249902959
125 x 175 x 240	6 boxes of 100 units	8720249902973
125 x 175 x 240	6 boxes of 100 units	8720249902997
125 x 175 x 240	6 boxes of 100 units	8720249903017
125 x 175 x 240	6 boxes of 100 units	8720249903031
125 x 175 x 240	6 boxes of 100 units	8720249903055



Medicom Healthcare B.V.  
Parallelweg 80a | 3931 MT Woudenberg | The Netherlands  
cs@medicom-eu.com | +31 332 866 460