

blackcurrant

aniseed

liquorice

mandarin

ART
by **RISKONTROL**[®]

Liven up your practice

RISKONTROL[®]

Infection control

- ART**
by **RISKONTROL**[®]
Non-sterile tips
-  Ref. 200292
 -  Ref. 200296
 -  Ref. 200297
 -  Ref. 200298
- CE
- Classic**
By **RISKONTROL**[®]
Non-sterile tips
-  Ref. 200280
 -  Ref. 200275
 -  Ref. 200290
 -  Ref. 200295
 -  Ref. 200285
- CE
- Sterile**
By **RISKONTROL**[®]
Sterile tips
-  Ref. 200255
 -  Ref. 200250
 -  Ref. 200265
- CE 0120
- Total Protec**
By **RISKONTROL**[®]
Non-sterile tips
with plastic sheath
-  Ref. 200300
 -  Ref. 200310
- CE

Professional dental products



*ART by RISKONTROL[®] is a registered trademark of PIERRE ROLLAND. All rights reserved. No information or part of this document may be reproduced or transmitted in any form without the prior permission of PIERRE ROLLAND.

NEW COLLECTION

RISKONTROL[®]

PIERRE ROLLAND
ACTEON Pharma

New: **4** colors, **4** fragrances

mandarin

aniseed


blackcurrant

liquorice

ART
by **RISKONTROL**[®]

Liven up your practice

Non-Sterile tips,
Box of 250 tips.



RISKONTROL®

Infection control

• Total protection

- protection for the patient and the practitioner,
- no danger of direct or cross-contamination,
- no danger of microbial migration in the unit's air and water circuits,
- disposable.

• Practical

- easy to fit and to remove,
- rounded tips with no risk of injury to the patient,
- flexible and bendable tips with non-shape-memory.

• Economical

- saves time (no need for cleaning/disinfection/sterilization, so, no need for traceability, etc).




• Ecological

- tips made from entirely recyclable food-quality plastic,
- environmentally-friendly,
- non-polluting in the case of combustion.

• Universal

- adapters available for any type of equipment,
- adapter in anodized aluminium or resin, very simple to screw or clip onto the syringe of your unit.



• Dry air

Due to its unique system, widely patented throughout the world, two independent tubes guide air and water from the syringe to the extremity of the tip, ensuring totally dry air.

Total Protec By RISKONTROL® For even more protection

Non-sterile tips with plastic sheath.
Box of 100 tips.
Available in white and blue.

Classic By RISKONTROL® The original

Non-sterile tips.
Box of 250 tips.
Available in white, blue, green, yellow and pink.

Sterile By RISKONTROL® The specialist

Sterile tips.
Box of 100 tips in individual blisters, sterilized with ethylene oxide.
Available in white and blue.

